USO PREVISTO:

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INTENDED USE:

Suitable for preventing adhesion between 3D resin-printed models and dental devices made of self-curing, heat-curing, and light-curing acrylic resins.

COMPOSITION:

The insulator composed of:

- Biodegradable Polymer
- Water
- Low molecular weight alcohols
- Radical inhibitor
- Plasticizer

DESCRIPTION:

Versatile use insulator, with a fast drying time, moderate viscosity and ideal for its application, has a translucent appearance to prevent contamination of dental structures. Specifically formulated to prevent chemical adhesion between 3D printed models and the structures to be fabricated on them, which are based on self-curing, heat-curing and light-curing acrylic resins, thus facilitating the separation process.

INSTRUCTIONS FOR USE:

For self-curing and light-curing structures

Apply a thin layer of insulator onto the 3D printed model using a soft-bristle brush, ensuring an even distribution across the entire surface. After application, allow the layer to dry for 15 to 20 minutes.

If a faster drying process is needed, gently apply compressed air for about 2 minutes, taking care not to deform or remove the release layer. Alternatively, a hair dryer on medium heat may be used for approximately 1 minute.

For heat-curing structures

Apply a thin layer of insulator onto the 3D-printed model using a soft-bristle brush, ensuring an even distribution over the entire surface. Allow it to dry for 5 minutes, then apply a second thin layer. Afterwards, let it dry for 15 to 20 minutes.

If faster drying is required, gently apply compressed air for approximately 2 minutes, taking care not to deform or remove the applied product. Alternatively, a hair dryer on medium heat may be used for about 1 minute.

Post-processing instructions:

Take care to remove the polymerized structure from the 3D printed model just after the polymerization process is finished, since if this is not done, there is a risk that the built structure will adhere to the 3D model.

To remove the fabricated structures on the printed model, it is suggested to use a rigid and thin spatula, with which a lever effort is generated to facilitate their separation.

In the specific case of self-curing resins, it is suggested that the spatula be introduced palatally/lingually so that the effort applied is more efficient, see image 1.

Image 1. Separation between the 3D model and a retainer made of self-curing acrylic.

WORKING IN CONJUNCTION WITH OTHER DEVICES:

This product is compatible with self-curing, heat-curing and light-curing acrylics. The instructions for use of each acrylic should be strictly followed to fabricate dental devices on 3D resin-printed models.

RESIDUAL RISKS:

-Inadequate storage can generate a change in the properties of the product, causing an increase in drying times, appearance of solid lumps and loss of product efficiency. This can be avoided with a proper storage of the product according to the mentioned indications.

-The improper application of this product can generate problems to insulate the structures manufactured on the 3D printed models, in addition to the application of the insulator using the finger or any other part of the body to spread it

on the structure of the model can cause allergies and irritation in the contact area.

-Inhalation of this product can cause dizziness, drowsiness and vertigo, so it is important to avoid leaving the product uncovered and exposed to the elements, as well as to avoid inhalation.

CONTRAINDICATIONS:

In the case of over-applying a very thick layer of the product, it may take longer to dry and additionally the layer formed on the 3D model will be thick enough to affect the morphology of the built acrylic structure. It is advisable to apply a thin and homogeneous layer on the entire surface of the model to avoid this inconvenience.

WARNINGS:

This product should not be used in an application other than that specified in the intended use. Do not inhale its vapors because you will experience dizziness and headache if you do.

PRECAUTIONS:

It should not be used by children, pregnant or breastfeeding women.

RECOMMENDATIONS:

Use a soft bristle brush to ensure that the applied layer of product is homogeneous and completely covers the surface of the model to be isolated.

STORAGE:

The product should be stored in its original container, at room temperature (< 30 °C) avoiding the following conditions:

- Direct exposure to sunlight or other types of radiation.
- High heat or humidity sources.
- Dust or any other type of contaminant.

FINAL DISPOSAL:

Do not dispose this product in water sources or sewage systems. This product should not be discarded in soil or subsoil. Disposal should be carried out through an authorized waste management company, in accordance with national and local regulations.

SHELF LIFE:

2 years.

	Consult instructions for use
	Maximum storage temperature
	Caution: Consult the accompanying documents. See the instructions for additional warnings or suggestions for proper use of the device.
	Keep away from sunlight
	Wear gloves
	Caution: The product may cause respiratory irritation, drowsiness, allergic skin reactions, and eye irritation
	Caution: Hazardous to the environment
	Keep in a dry place
	Eye protection

 Fabricado por/Manufactured by
New Stetic S.A. Cra 53 N° 50-09 Guarne - Antioquia-Colombia.
Tel: (+ 57-604) 403 8760 - 550 0000
infocolombia@newstetic.com / www.newstetic.com