USO PREVISTO:

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INTENDED USE:

PORTUX 3D TRY-IN is a material specifically formulated for the 3D printing of full or partial dental structures intended for short-term clinical try-in procedures. It is designed to accurately assess design fit, aesthetics, phonetics, and occlusion during trial appointments. This material must be handled exclusively by dental professionals or dental technicians who are properly trained in the use of photopolymerizable resins for 3D printing.

- Crown and Bridge Restorations specific for design evaluation
- Monolithic Denture Prosthetics specific for design evaluation
- Screw Retain Denture Prosthetics specific for design evaluation
- Hybrid Denture Prosthetics specific for design evaluation

COMPOSITION:

The resin is a mixture of:

- Acrylic monomers.
- Photoinitiators (Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide).
- Stabilizing additives.
- Pigments and dyes.

DESCRIPTION:

Low-viscosity, photopolymerizable resin available in tooth shades, suitable for the fabrication of provisional full or partial dental structures for try-in purposes. It offers excellent dimensional accuracy and high definition, faithfully reproducing the digital design. The product is compatible with open-system LCD/DLP stereolithography 3D printers operating at wavelengths between 385 and 405 nm.

INSTRUCTIONS FOR USE:

Mix the product for at least 20 minutes using a mechanical shaker or roller before first opening. This ensures proper printer performance and consistent color reproduction.

To minimize bubble formation, agitate the product at least 1 hour prior to use.

Print with the resin following your printer's handling and operating instructions.

Post-Processing instructions:

- Clean printed parts in isopropyl or ethyl alcohol (>90%), using an ultrasonic bath to facilitate cleaning. Immerse the prints in a container with alcohol for 5 minutes. Then immerse them in fresh alcohol for the same amount of time. It is recommended to apply compressed air gently between cleanings to remove excess resin from cavities or critical areas of the print.
- Remove the pieces from the alcohol and dry them, either by gently applying compressed air or in an oven at 40°C for 30 minutes. **IMPORTANT:** Avoid curing wet or damp prints, as this affects the accuracy and final definition of the printed part.
- Nitrogen post-curing units inhibit secondary chemical reactions, thereby ensuring enhanced color stability of the product. However, conventional units are adequate for curing PORTUX 3D Try-in resins (in accordance with the manufacturer's instructions) during short patient try-in periods.
- If using photopolymerizable stains and/or glazes, the curing times specified by the manufacturer must be applied during either the first or second exposure, without adding to the previously defined exposure times.
- Consult the post-curing times for the different equipment with the technical support team for PORTUX 3D resins or on the New Stetic USA website.

Note: Exposing the printed structure to curing times longer than specified may increase the material's fragility and promote the occurrence of fractures.

WORKING IN CONJUNCTION WITH OTHER DEVICES:

This product is compatible with photopolymerizable resin stains and glazes. For their application on structures printed with this resin, it is essential to follow the instructions for use provided by the product manufacturer.

RESIDUAL RISKS:

- Contact with uncured material or the release of residual unpolymerized components may cause allergic reactions. This risk is minimized by performing post-curing under the specified time and conditions. Biocompatibility testing ensures the safety of the device during its intended use.
- Printed try-in structures can be characterized with photopolymerizable resin stains or glazes. Contact with these materials before complete polymerization may trigger allergic reactions. To prevent this, polymerization must be

carried out according to each manufacturer's instructions.

- During use, printed dental structures are subjected to mechanical forces that may result in fractures. This risk is reduced by implementing an adequate design in the corresponding CAD software and by properly handling the printed object.
- Improper storage may alter the resin's polymerization properties. Exposure to light or other radiation sources can cause premature polymerization, rendering the material unusable. To avoid this, store the product according to the manufacturer's recommendations and always return it to its original container immediately after use.
- Misusing the printed pieces as provisional or definitive structures may lead to premature fractures, allergic reactions, or other adverse effects in the patient due to overexposure to the material. This risk is eliminated by using the try-in structures in the patient's oral cavity only for the stipulated time.

CONTRAINDICATIONS:
The use of dental test frameworks made of this material is contraindicated in patients with proven allergic reactions to products derived from methacrylate resins.

WARNINGS:
This product should only be used in the applications specified for its intended use. PORTUX 3D TRY-IN printed test structures are only biocompatible following the post-processing instructions.

***Not to be used as a temporary structure.**

PRECAUTIONS:
Use in children, pregnant or lactating women is not recommended.

STORAGE:
The product should always be kept within its original packaging, at room temperature (< 30°C) avoiding the following conditions:

- Direct exposure to sunlight or other types of radiation.
- High heat or humidity sources.
- Dust or any other type of contaminant.

FINAL DISPOSAL:
Do not dispose of the product into bodies of water, sewage systems, soil, or subsoil. Disposal must be carried out through an authorized waste management company, in compliance with applicable national and local regulations.

SHELF LIFE:
2 years.

	Consult instructions for use
	Maximum storage temperature
	Caution: Consult the accompanying documents. See the instructions for additional warnings or suggestions for proper use of the device.
	Keep away from sunlight
	Wear gloves
	Caution: The product may cause respiratory irritation, drowsiness, allergic skin reactions, and eye irritation
	Caution: Hazardous to the environment
	Keep in a dry place
	Eye protection

 Fabricado por/Manufactured by
New Stetic S.A. Cra 53 N° 50-09 Guarne - Antioquia-Colombia.
Tel: (+ 57-604) 403 8760 - 550 0000
infocolombia@newstetic.com