

Dientes de resina acrílica
con fluorescencia y cross link
Dos capas cromáticas

[Cross linked fluorescent
acrylic resin teeth
Two chromatic layers]


DIENTES DE RESINA ACRÍLICA

Olympic[®]


ACRYLIC RESIN TEETH



ANTERIORES SUPERIORES Upper ANTERIORS



ANTERIORES INFERIORES Lower ANTERIORS



3M
222
1N
2N
2P
4P
267
3R
4N

layers
[capas]

		3M A 9.0 B 7.0 C 38.5
		222 A 7.6 B 6.8 C 35.6
		1N A 9.8 B 8.3 C 44.5
		2N A 10.3 B 7.5 C 42.0
		2P A 11.3 B 8.5 C 47.0
		4P A 11.5 B 8.4 C 45.3
		267 A 12.4 B 9.5 C 51.0
		3R A 12.7 B 9.2 C 51.0
		4N A 9.4 B 8.0 C 42.2

		3M A 8.2 B 5.0 C 31.5
		222 A 8.2 B 5.0 C 31.5
		1N A 9.3 B 5.0 C 34.0
		2N A 9.0 B 5.0 C 33.0
		2P A 9.8 B 5.5 C 36.0
		4P A 11.5 B 5.5 C 38.0
		267 A 11.5 B 5.5 C 38.0
		3R A 11.5 B 5.7 C 39.2
		4N A 9.6 B 5.2 C 34.4

2D
3D
3N
1H
136
137

		2D A 9.5 B 7.8 C 42.0
		3D A 9.3 B 8.0 C 44.0

		2D A 9.2 B 5.0 C 32.0
		3D A 8.8 B 5.0 C 34.0



3N
A 10.0
B 8.0
C 45.0

1H
A 12.0
B 8.5
C 48.0

I36
A 10.2
B 9.0
C 48.4

I37
A 10.0
B 9.2
C 49.8



3N
A 9.5
B 5.0
C 34.0

1H
A 10.8
B 5.8
C 38.5

I36
A 9.3
B 5.5
C 37.5

I37
A 10.8
B 5.8
C 38.5



A23-223
A24-224
A25
3P
4H
264
A26



A23-223
A 8.1
B 7.5
C 38.5

A24-224
A 8.2
B 8.0
C 39.8

A25
A 9.5
B 8.5
C 46.0

3P
A 10.5
B 8.5
C 48.0

4H
A 9.4
B 8.2
C 44.3

264
A 10.8
B 8.0
C 45.0

A26
A 9.8
B 9.0
C 47.0



A23-223
A 8.8
B 5.0
C 34.0

A24-224
A 8.8
B 5.0
C 34.0

A25
A 9.5
B 5.5
C 34.0



3P
A 10.0
B 6.0
C 38.0

4H
A 9.3
B 5.5
C 37.5

264
A 10.0
B 6.0
C 38.0

A26
A 8.5
B 5.5
C 35.3

POSTERIORES
SUPERIORES E INFERIORES
UPPER-LOWER POSTERIOR



29Z
0°

A 27.3
C 6.5
B 28.8
D 7.0



29L
20°

A 29.0
C 8.2
B 32.5
D 8.2



31Z
0°

A 29.5
C 7.2
B 31.0
D 8.4



31S
20°

A 30.0
C 8.2
B 34.5
D 8.2



32F
10°

A 31.0
C 8.5
B 32.5
D 9.0



31M
20°

A 31.0
C 7.9
B 32.5
D 8.5



33F
10°

A 32.0
C 9.1
B 34.5
D 9.4



31L
20°

A 30.9
C 8.7
B 34.5
D 8.7



29S
20°

A 29.0
C 7.5
B 31.5
D 8.2



33M
20°

A 32.0
C 8.6
B 35.5
D 8.7



29M
20°

A 29.0
C 7.2
B 32.0
D 7.8




33L
20°

A 33.5
C 9.3
B 36.5
D 9.3



ACRYLIC RESIN TEETH[®] Olympic[®]



233L

20°

A 32.0
C 10.1
B 32.2
D 9.8



PBS

20°

A 28.7
C 8.1
B 31.4
D 8.2



PDS

20°

A 30.2
C 8.3
B 32.6
D 8.5



PF

20°

A 31.2
C 8.9
B 36.1
D 9.4



30M

33°

A 29.0
C 7.7
B 31.5
D 7.7



30L

33°

A 28.7
C 8.5
B 31.1
D 8.2



32M

33°

A 31.0
C 7.6
B 33.5
D 8.2



32L

33°

A 29.4
C 9.3
B 32.7
D 9.5



34M

33°

A 32.5
C 8.8
B 34.5
D 9.3



34L

33°

A 32.6
C 9.8
B 32.5
D 10.0



230L

33°

A 29.5
C 9.0
B 29.8
D 8.8



Escanea estos códigos para/scan these codes to:

Visitar nuestra página en facebook
Visit our page on facebook




Descargar esta carta de moldes en PDF
Download this mould chart in PDF



Articulación [ARTICULATION]

ANTERIORES [ANTERIORS]	0°	10°	20°	33°
3M	29Z		29S-29M	30M-30L
2D	31Z		233L-31M	32M-32L
3D			31M-31S	32M-32L
222			29S-29M	
223	29Z		29L-29M	30M-30L
A23	29Z		29L-29M	30M-30L
A24			29S-29M	30M
A25	29Z		31L-PBS	230L
A26		32F	31S-PDS	
1N			31S-31L-PBS	
2N		32F	233L-31M	32L
3N	31Z		29L	32L-230L
4N			31S	32M-32L
2P			31L-PDS	32L
3P	31Z		PDS	32L
4P			31S-31M	32L
1H			233L-33L-33M	34M-34L
4H		32F	233L-31M	230L
136	31Z		233L-31M	230L
137		33F	233L-PF	34M
264	31Z		PDS	32L
267			33L-33M	34M-34L
3R		33F	33L	34M-34L


Angulación de posteriores [POSTERIOR ANGLATION]



Colores [SHADES]

03 • 59 • 62 • 65 • 66 • 67 • 69 • 77 • 81


Ejemplos [EXAMPLES]



Mould 3M

ANTERIORES SUPERIORES [UPPER ANTERIORS]

A 9.0
B 7.0
C 38.5



Mould 3M


ANTERIORES INFERIORES [LOWER ANTERIORS]

A 8.2
B 5.0
C 31.5

A LARGO DEL CENTRAL SIN CUELLO (mm)
[LENGTH OF CENTRAL WITHOUT COLLAR]

B ANCHO DEL CENTRAL (mm)
[WIDTH OF CENTRAL]


C ANCHO DE 6 ANTERIORES EN LÍNEA RECTA (mm)
[WIDTH OF 6 ANTERIORS ON A STRAIGHT LINE]



Mould 29S

POSTERIORES SUPERIORES [UPPER POSTERIORS]

A 29.0
C 7.5



Mould 29S

POSTERIORES INFERIORES [LOWER POSTERIORS]

B 31.5
D 8.2

A-B ANCHO DE 4 POSTERIORES EN LÍNEA RECTA (mm)
Width of 4 posteriors on a straight line

C-D ANCHO DEL PRIMER MOLAR SUPERIOR O INFERIOR IZQUIERDO
Width of the upper first molar or lower left molar

ACRYLIC RESIN TEETH
Olympic®

Síguenos / Follow us:



New Stetic S.A.
Cr. 53 No.50-09 Guarne, Antioquia, Colombia
PBX: (+57-604) 550 00 00 - 403 87 60
infocolombia@newstetic.com
www.newstetic.com

ISO
22112



SC 358-1
SC 358-2



DM 017-1
DM 017-2



REV. 22-06-2023