


Reflection[®] Acrylic teeth Dientes acrílicos

Carta de moldes
MOULD CHART



Estética Natural
Natural Esthetics

Resistencia a la Abrasión
Abrasion Resistance

Four
chromatic **layers**
Cuatro capas cromáticas




Dientes acrílicos

Reflection®

Acrylic teeth

Anteriores Superiores / Upper anteriors




RO25

A 9.6
B 8.0
C 41.7




RT36

A 9.9
B 7.9
C 42.9




RO98

A 11.3
B 8.9
C 47.5




RT56

A 9.7
B 8.2
C 43.9



RO99

A 12.5
B 9.7
C 51.4




RT66

A 10.0
B 8.0
C 42.9

4 Four
chromatic layers

Cuatro capas cromáticas



RT76

A 10.2
B 8.1
C 44.8




RT77

A 10.3
B 8.5
C 45.3



RT88

A 10.5
B 8.6
C 46.1



Estética Natural

Natural Esthetics

Anteriores Inferiores

Lower anteriors



R3
A 8.0
B 5.1
C 33.1



R5
A 9.3
B 5.1
C 34.4





R7
A 8.4
B 5.4
C 35.3



R13
A 8.6
B 5.6
C 36.2



R14
A 10.6
B 6.1
C 41.4



Posteriores superiores-inferiores

Upper-Lower posteriors



RC41
X 29.2
U 8.8
Y 31.0
L 8.3




RC45
X 34.5
U 9.4
Y 35.9
L 8.4



RS11
X 29.5
U 8.1
Y 31.3
L 7.1



RS13
X 30.4
U 8.5
Y 33.4
L 7.3



ARTICULACIÓN


Articulation

Anteriores / Anteriors
Posteriores / Posteriors


RO25	R3	RC41-RS11
RO98	R14	RC45
RO99	R14	RC45
RT36	R7	RC41-RS11
RT56	R5	RC41
RT66	R7	RC41
RT76	R5	RC41-RS11
RT77	R5	RS13
RT88	R13	RS13

Ejemplos / Examples

Molde / Mould RT36
[Anteriores Superiores]
[Upper Anteriors]




Molde / Mould R7
[Anteriores Inferiores]
[Lower Anteriors]




- A** Largo del central sin cuello (mm)
Length of central without collar (mm)
- B** Ancho del central (mm)
Width of central (mm)
- C** Ancho de 6 anteriores en línea recta (mm)
Width of 6 anteriors on flat plane (mm)

Molde / Mould RS13
[Posteriores Superiores]
[Upper Posteriors]



Molde / Mould RS13
[Posteriores Inferiores]
[Lower Posteriors]



- X-Y** Promedio del ancho de un juego 1x4
The average width of a 1x4 set
- Y-L** Promedio del ancho del primer molar superior o primer molar izquierdo
The average depth of the upper first molar or lower left molar

Colores / Available Shades:

A0 - A1 - A2 - A3 - A3.5 - A4 - B0 - B1 - B2 - B3 - B4 - C1 - C2 - C3 - C4 - D2 - D3 - D4



Síguenos / Follow us:



ISO
22112



SC 358-1
SC 358-2



DM 017-1
DM 017-2



New Stetic S.A.
Cr. 53 No.50-09 Guarne, Antioquia, Colombia
PBX: (+57-60-4) 550 00 00 - 403 87 60
FAX: (+57-60-4) 551 31 34
infocolombia@newstetic.com - www.newstetic.com